Print Management System Printlinkの開発

The Development of the Printlink Print Management System

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Yamamichi, Yoji

Printlink System is a hospital medical image transfer system binding between medical image diagnostic devices such as CT, MRI or Regius and Laser imagers or image archiving systems in the network environment.

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
Fig. 2 Block diagram of Printlink System

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		Image Capture for DICOM Network Printlink-IV / ID / DV / HG / IN	DICOM Print Server Printlink-OD
Image Data Input CH		Max. 1ch	Max. 16ch
Image Data Output CH		Max. 8ch (Selectable)	Max. 1ch
Interface	Input	Video / Digital / Digital Video High Speed Video / Ethernet (DICOM, Toshiba Ethernet Protocol)	Ethernet 10/100BaseTX supports DICOM Basic Print SCP extended for Konica Imager functions.
	Output	Ethernet: 10/100BaseTX supports DICOM Basic Print SCU extends for Konica Imager functions. And DICOM Storage SCU	RS485 for Data Transfer and RS-232C for Imager control
ID Input		OCR or Touch panel Keypad	
Image Depth		8 bit or 12 bit	8 bit or 12 bit
Imager		DICOM Printer Li-1417D/HDH	Li-10A, Li-21, Li-77A/8, Li-62P
Image Spool Memory		4GB	4GB
Host Control		Available using major camera control Protocol	
Print functions		All print functions in Konica's Imager can be set by the Keypad of Printlink.	All functions from Printlink can be supported. Depending on the supported function in case of DICOM SCU.
Operational Keypad		Touch panel on the LCD with Back Light	

- Printlink Input Video with pixel clock from 3.75-60MHz
- Printlink ID Input through RS485 in 8 or 12bit
- Printlink HG Input Video with pixel clock from 80-135MHz
- Printlink DV Input Digital Video
- Printlink IN Input through Ethernet, DISCOM, Toshiba Protocol
- Printlink OD Output Digital for Lazer Imager

Fig. 3 KONICA Printlink Specifications

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