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Automatic Detection of DR Images in which the Lungs are Partially Indistinguishable
and of DR Images that are Motion Blurred

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Abstract

Among hospital radiologic technologists who must make the rounds of patients unable to go to the x-ray room, x-ray photography using DR (digital radiography) has become popular. Konica Minolta's portable AeroDR solution realizes work flow efficiency because the AeroDR provides immediate image access: images can be viewed in the same room where the images of the patient are taken. Thus, if an image taken proves inappropriate for diagnosis, another image can be taken immediately.

However, this potential depends upon the ability of the radiologic technologist to properly judge the quality of images taken. In front chest images, which are frequently taken at the patient's bedside, it can be difficult to keep the patient motionless during image capture or to confirm proper cassette positioning for an x-ray exposure area. Therefore, there are cases of DR images in which the lungs are partially indistinguishable as well as DR images that are motion blurred. Since it is difficult to differentiate such images from appropriate ones, the radiologic technologist may fail to retake an x-ray when needed.

With the goal of reducing the risk of overlooking diagnostically inappropriate chest images, we developed functions to automatically detect images in which the lungs are partially indistinguishable and images that are motion blurred. We implemented these functions in the portable CS-7 console, which functions alongside the AeroDR. Identifying images in which the lungs are partially indistinguishable is accomplished by extracting effective features that differentiate appropriate images from inappropriate images, and then by judging images taken with a classifier that employs a vast degree of knowledge derived from a huge database. With motion blurred images, detection is based on features of both motion blurred and static images, which are extracted after analyzing signal components of lung fields between the ribs. Presented here are these detection functions and the technologies used therein.

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
Fig. 1 Example of a DR image in which the lungs are partially indistinguishable. The majority of times when x-rays must be retaken are instances when images in which the lungs are partially indistinguishable due to a slight shift of patient position during image capture. In this example, if there is an existent tumor mass, failure to distinguish the lungs results in failure to observe that tumor mass (see arrow).


Fig. 2 Example of a motion blurred DR image. The blurred portion (right image) is often missed in a whole lung image (left image) because of the low resolution of the CS-7 console's screen.

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
Fig. 3 Example of a CS-7 console screen image in which partially indistinguishable lungs (red rectangle) and motion blurring (orange square) are automatically detected.

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
Fig. 4 Flow chart of an image detection algorithm for the detection of partially indistinguishable lungs. Effective features are extracted from various areas of an image to differentiate an image in which the lungs are partially indistinguishable from a normal image. Whether or not the lungs are partially indistinguishable is determined by a classifier that employs a vast degree of knowledge derived from a huge database.

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
Fig. 5 Flow chart of an image detection algorithm for the detection of motion blurred images. A plurality of small areas are set in a lung field between the ribs, and features of both a motion blurred image and a static image are extracted by analyzing signal components of each small area. The detection of motion blurred images is achieved based on these extracted features.

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