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Development of Bone Suppression Processing as an Application of
Computer Aided Detection of Nodules in Chest Radiographs

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Abstract

Chest x-ray images are widely used for their convenience, but interpreting chest images is difficult. Images of many components such as bones, blood vessels, and the heart are superimposed, and a wide variety of conditions must be diagnosed. Because of this, chest images are a field in which computer-aided diagnosis is particularly beneficial. To assist interpretation, we developed bone suppression processing that suppresses signals of ribs and clavicles in chest images. This lessens the chance of overlooking conditions due to subtle chest images, and it reduces false positives.

Twenty doctors participated in an observer performance study to evaluate the effectiveness of our bone suppression processing. In the study, we compared performances in 40 normal cases and 40 nodule cases (≤ 30 mm). Using bone suppression processing, the average area under the receiver operating characteristic (ROC) curve improved from 0.87 to 0.90 with statistical significance. We grouped the doctors into two groups based on diagnosis experience. Doctors with ten or fewer years experience benefitted the most from bone suppression processing.

In this paper, we give an overview of our bone suppression algorithm, which can robustly and selectively suppress signals of ribs and clavicles. We also present our new system, in which original and bone suppressed images can be viewed interchangeably without changing workflow. We also introduce examples of incorporating the bone suppression system with our console system, Unitea α , and we discuss our PACS viewer, I-PACS EX.

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**Konica Minolta Laboratory U.S.A., INC. Health Care Division

Table 1 Dual energy subtraction and bone suppression.

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	Dual energy subtraction (Two shot method)	Bone suppression
Purpose	To improve visibility of lesions behind bones and bone lesions.	To improve visibility of lesions behind ribs and clavicles.
Method	Bone and soft-tissue images are made by two exposures, one of low voltage and one of high voltage, based on the patient's permeability of radiation.	Ribs and clavicles are recognized from their anatomical features using image recognition technology, allowing signals from the bones to be suppressed.
Price	Expensive because a special exposure device is required.	Inexpensive because no special exposure device is required.
Earlier image	Inapplicable to earlier images.	Applicable to earlier images.
Exposure	Large amounts of exposure due to two exposures.	Equivalent to a single regular exposure.
Image quality	Due to two exposures, motion blur occurs. Image quality is susceptible to body size of the patient.	Bone recognition errors are seldom. Normal tissue running in the same direction as accompanying bone may possibly be suppressed.

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Original image

Bone suppressed image

Fig. 1 A posterior-anterior chest image with a nodule (above) and a bone suppression image of the posterior-anterior chest image (below).

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
Fig. 2 The flow of bone suppression processing.

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
Original image

Bone suppressed image

Fig. 3 A nodule overlapped with a posterior rib.

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Table 2 Observer detection study.

Item	Content																
Method	<ul style="list-style-type: none"> Sequence test with a continuous confidence-rating scale Viewer: ROCViewer V11.4.0.7¹⁰⁾ 																
Analysis	<ul style="list-style-type: none"> ROC (receiver operating characteristic) analysis Rockit V0.9.1 & DBM MRMC V2.2 LROC (localized ROC) analysis 																
Observers	<table border="1"> <thead> <tr> <th>Class</th> <th>Surgery</th> <th>Internal medicine</th> <th>Total</th> </tr> </thead> <tbody> <tr> <td>Pulmonology residents</td> <td>6</td> <td>7</td> <td>13</td> </tr> <tr> <td>Experienced pulmonologists</td> <td>3</td> <td>4</td> <td>7</td> </tr> <tr> <td>Total</td> <td>9</td> <td>11</td> <td>20</td> </tr> </tbody> </table>	Class	Surgery	Internal medicine	Total	Pulmonology residents	6	7	13	Experienced pulmonologists	3	4	7	Total	9	11	20
	Class	Surgery	Internal medicine	Total													
	Pulmonology residents	6	7	13													
	Experienced pulmonologists	3	4	7													
Total	9	11	20														

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
Fig. 4 ROC curves for detection of nodules by all observers. Detection capability is higher with bone suppression processing than without.

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
Fig. 5 Receiver operating characteristic curves for detection of nodules by experienced pulmonologists and inexperienced residents. With each, detection is higher with bone suppression processing.

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
Fig. 6 Comparison of area under the curve for the localized ROC curves by all observers. In the interpretation of radiographs using bone suppression processing, doctors with less diagnosis experience (B) improved more than those with more experience (A).

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Fig. 7 The Unitea α console system, a product in which bone suppression processing is incorporated. The Unitea α is Konica Minolta's original all-in-one solution, combining reliability and superior throughput.

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
Fig. 8 Because the original image and the bone suppressed image are captured from the same viewpoint, the radiologist can readily discern what bone suppression processing reveals.

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Fig. 9 BS processing will have been used with I-PACS EX, in which a rapid diagnostic environment is realized with various monitors that provides abundant functionality and simple operation.

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