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Development of a Prostate Cancer Diagnostic System Based on Surface Plasmon Field-Enhanced Fluorescence Spectroscopy

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Abstract

With the ageing of society and the Westernization of dietary habits, incidences of prostate cancer (PC) have significantly increased in the past few years. PC begins in a man's prostate gland. Prostate specific antigen (PSA) screening helps identify early-stage cancer, when treatment is most effective. However, PSA lacks cancer specificity, as it can also be elevated in conditions such as benign prostatic hyperplasia (BPH) and chronic prostatitis. Due to the low specificity of PSA screening and the resultant high number of false positives, a large number of unnecessary biopsies are performed each year in order to obtain definite diagnoses.

Just such a cancer-specific marker may be the structural change in carbohydrate chains produced by cell changes during carcinogenesis. That same cancer-specific change was also found on PSA, which was conventionally detected by using column chromatography combined with the ELISA method. This satisfies all specificity requirements for quantitative analysis of a PSA-containing, cancer-specific carbohydrate (LacdiNAc-PSA) in serum. However, it is impractical for clinical use because it takes an excessive amount of time and labor.

We have developed an automated immunoassay system employing surface plasmon field-enhanced fluorescence spectroscopy (SPFS), one of the most sensitive fluorescence detection technologies today. The performance of LacdiNAc-PSA immunoassay with SPFS showed high sensitivity. Clinical evaluation revealed that the level of LacdiNAc-PSA in serum was significantly higher in cases of PC than in cases of BPH. Further, the level of LacdiNAc-PSA in serum correlated positively with the progression of cancer. Thus, a PSA carbohydrate marker in serum may eliminate the need for a needle biopsy in cases of suspected prostate cancer.

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
Fig.1 The principle of surface plasmon field-enhanced fluorescence spectroscopy (SPFS).

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
Fig. 2 The automatic SPFS immunoassay system with a disposable plastic sensor chip.

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
Fig. 3 (a) Schematic view of sandwich immunoassay on an SPFS chip and ELISA plate, and (b) calibration curves of total PSA in female serum derived from SPFS (●) and ELISA (Δ). The SPFS immunoassay has both higher sensitivity and a wider measurement range than those of ELISA.

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Table 1 Performance of six automated immunoassay systems that perform total PSA tests.

Instrument	Company	Assay method	Measurement range (pg/ml)
—	Konica Minolta	SPFS	0.04 – 100000
ARCHITECT	Abbott	CLIA	8 – 100000 ⁵⁾
Unicel	Beckman	CLIA	8 – 150000 ⁶⁾
IMMULYZE	Siemens	CLEIA	5 – 20000 ⁷⁾
ECLusys	Roche	ECLIA	2 – 100000 ⁸⁾
AIA	Tosoh	FEIA	10 – 100000 ⁹⁾

CLIA: chemiluminescence immunoassay
 CLEIA: chemiluminescence enzyme immunoassay
 ECLIA: electro chemiluminescence immunoassay
 FEIA: fluorescence enzyme immunoassay

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
Fig. 4 (a) Lectin sandwich immunoassay with SPFS and (b) calibration curves of LacdiNAc-PSA using lectin sandwich immunoassay with SPFS (●) and with ELISA (△). The limit of detection (LOD) of LacdiNAc-PSA using automated SPFS was 20.0 pg/ml. The LOD of LacdiNAc-PSA using ELISA was 20 ng/ml. The automated SPFS immunoassay clearly had higher sensitivity than that of the ELISA immunoassay.

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
Fig. 5 (a) Measured values of total PSA using ELISA and (b) plot of the measured values of LacdiNAc-PSA using lectin sandwich automated SPFS immunoassay for prostate cancer (PC) (N = 39) and benign prostate hyperplasia (BPH) (N = 31) in the gray zone (total PSA concentration: 4 – 20 ng/ml).

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
Fig. 6 Receiver operating characteristic (ROC) curve analysis for LacdiNAc-PSA (●) and total PSA (●). The area under the ROC curve for LacdiNAc-PSA (0.823) was significantly greater than the area under the total PSA curve (0.506).

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Table 2 Diagnostic serum marker performance. In the 4 to 10 ng/ml PSA range at 90% and 95% sensitivity, both the specificity and the area under the ROC curve (AUC) of the LacdiNAc-PSA exceeded those of total PSA, %fPSA, and Phi. (AUCs LacdiNAc-PSA 0.866, PSA 0.516, %fPSA 0.648 and Phi 0.708).

Markers	95% Sensitivity		90% Sensitivity		ROC analysis AUC
	Cut off	%Specificity	Cut off	%Specificity	
LacdiNAc-PSA	27650	52.0	32800	66.0	0.866
total PSA	4.2	7.5	4.4	10.8	0.516
%fPSA	27.9	9.9	24.1	19.8	0.648
Phi	22.1	14.1	27.0	31.1	0.708

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
Fig. 7 LacdiNAc-PSA serum levels correlated with tumor Gleason scores (GS) of PC patients. The LacdiNAc-PSA serum levels were significantly higher in PC (GS \geq 7) serum than in PC (GS \leq 6) serum (P=0.0350): (a) Plot of the measured values of total-PSA using ELISA and (b) plot of the measured values of LacdiNAc-PSA using SPFS immunoassay for PC (GS \geq 7 (N=14), GS \leq 6 (N=10)), BPH (N=14), and prostatitis (N=7) in the gray zone (total PSA concentration: 4–20 ng/ml), healthy men control (HMC) (N=4).

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