# HSTT (High Sensitive Tissue Testing) 技術の紹介

HSTT (High Sensitive Tissue Testing) Technology

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
#### **Abstract**

Konica Minolta succeeded in developing a digital pathology technology HSTT at molecular imaging levels. The HSTT technology is based on Konica Minolta's heritage of photographic chemistry. We are accumulating data which indicate that the HSTT technology can detect antibodies which cannot be detected by an immunohistochemistry (IHC) such as the DAB method.

The HSTT is a breakthrough in digital pathology which uses immunohistochemistry to enable both antigen quantitation and localization (visualization). Many researchers have attempted to apply immunofluorescence to protein quantitation, but they have faced difficulties in attaining sufficient brightness and stability. Konica Minolta successfully developed innovative nanoparticles whose brightness exceeds the auto-fluorescent noise level in tissue samples and which have excellent solvent durability during the pathology staining process. Furthermore, with our image recognition and digital processing technologies, we have achieved nanometer-level resolution in quantifying the number of particles (protein expressions).

The HSTT technology is characterized by (1) a high sensitivity, (2) a quantitative and wide dynamic range, (3) a simultaneous detection of multiple types of proteins, and (4) an analysis of locations of detected proteins in cells and tissues. In principle, since the nanoparticles can bind to any antibodies through surface modification, we succeeded in staining dozens of antibodies. The HSTT technology can also be applied to detection of mRNA, miRNA, and DNA, cytology, and blood testing, and we are developing research collaborations in these areas as well.

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

Fig. 1 Basic properties of PIDs (phosphor integrated dots).

(Fig. 1a) A scanning electron microscopy image of PID nanoparticles. (Fig. 1b) Comparison of 580 nm light-excited fluorescence intensities among Dyes (sulforhodamine 101), QDs (quantum dots: Qdot 625), and PIDs. Each fluorescence intensity is shown on a logarithmic scale on the vertical axis. (Fig. 1c) Comparison of photostabilities of Dyes and PIDs. The Dyes or PIDs on a glass slide were irradiated with 580 nm light, and the fluorescence intensity of both fluorescent materials was measured over time. The fluorescence intensity of Dyes rapidly decreased to 2% of the initial value in 10 min. whereas that of PIDs remained over 80% of the

initial value under the same irradiation condition.

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



Fig. 2 Dynamic range of PIDs imaging.

Sensitivity of imaging of PIDs (Fig. 1d) was compared with that of DAB (Fig. 1a), QDs (Fig. 1b) and Dyes (Texas Red) (Fig. 1c). Biotin solutions of various concentrations were incubated with deparaffinized paraffin sections of human breast cancer tissues to prepare mimic slide samples. After washing them, streptavidin-conjugated HRPs visualized by DAB, streptavidin-conjugated Dyes (Texas Red), streptavidin-conjugated QDs, and streptavidin-conjugated PIDs were reacted with the biotins which were bonded on each slide glass. Streptavidin-conjugated DAB was analyzed using Aperio image analysis system (Leica) to quantify the degree of brown color of IHC-DAB. DAB showed insufficient quantitative capability. Each of QDs and Dyes showed a narrow dynamic range due to weak fluorescent intensity, whereas PIDs showed a high sensitivity and a linear wide dynamic range between 10<sup>-5</sup> to 10<sup>-3</sup> mM of biotin concentration because of strong brightness of each particle. (HRPs: horseradish peroxides)

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
Fig. 3 Quantitative performance of PID imaging.

Cultured cells of eight types of breast cancer which expressed HER2 protein at various levels were selected to investigate whether PIDs can accurately measure the protein expression level in cultured cancer cells. These cells were immunostained with monoclonal anti-HER2 antibody, biotinylated monoclonal secondary antibody, and PIDs. The HER2 expression levels of these cells were also examined using ELISA which is a conventional method of quantifying protein. The figure shows a very high correlation.

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
Fig. 4 A double stained image of ki67 protein (green) and HER2 protein (red) in cultured cells.

As shown in this photo, the technology using the PID nanoparticles can detect multiple markers on an identical tissue by altering fluorescent colors.

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