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Clinical Value Proposal of Ultrasound Diagnostic System SONIMAGE HS1

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Abstract

Ultrasound diagnostic equipment provides non-invasive, real-time medical imaging and is widely used in clinical practice. In 2014, Konica Minolta launched its first ultrasound system, the SONIMAGE HS1, whose features included 1) expanded areas of the body that can be imaged, 2) simple and customizable operation, and 3) quick and reliable network connection. Users rated these features highly, and user feedback guided us in further refining these features. Thus, in 2016, we created even more advanced technologies in the latest version of the SONIMAGE HS1.

To expand the areas of the body that the SONIMAGE HS1 can image, we developed a new, 7 MHz wideband linear array probe, the L11-3. This probe accesses a wide range of depth in musculoskeletal and superficial applications, including targets located less than 1 cm beneath the surface of the skin. Additionally, the probe employs a new method of detecting structures as tiny as several hundred microns.

The range and simplicity of operating the SONIMAGE HS1 are supported by specialized probes. In addition to the L11-3, the newly developed HL18-4 hockey-stick linear probe deals with complex procedures and puncturing, while the MC10-3 micro-convex probe allows examination of deeper areas of the body from a narrow gap. In addition, the monitor's screen layout can be customized to a user's unique workflow.

The clean and easy network connectivity of the SONIMAGE HS1 allows patient ID information to be automatically shared with an electronic health record (EHR) system at the click of a button, with no need to manually key in ID data.

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
Fig. 2 Linear probe structure.

The multi-level acoustic matching layer matches the acoustic impedances of the piezoelectric element with that of a living organism.

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
Fig. 3 Simulated frequency spectrum for various impedances of the first level of the acoustic matching layer.

The bottommost level of the of the acoustic matching layer borders the piezoelectric element. This level, designed via computer simulation, optimizes impedance matching between the piezoelectric element and the acoustic matching layer.

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
Fig. 4 The -6dB bandwidth obtained in Fig. 3 and the acoustic impedance of the bottommost acoustic matching layer level.

The -6dB bandwidth is almost saturated at 15 MRayls of acoustic impedance. The impedance of the first acoustic matching layer level is determined on the basis of this simulation result.

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Fig. 5 Wideband linear probe L11-3.

The L11-3 probe has a view field of 40 mm and accesses from shallow to deep tissues.

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
Fig. 6 Clinical L11-3 probe images: (a) Femoral nerve block procedure with a depth setting of 2.5 cm. (b) Sciatic nerve block procedure with a depth setting of 6.0 cm. The L11-3 probe accesses both visual field depths.

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
Fig. 7 Frequency characteristics of signals received by the piezoelectric transducer.

In conventional signal processing (a), the signal received is not filtered and has a Gaussian characteristic. In our new signal processing (b), the signal received is filtered to attain a flat frequency response. This flat frequency response provides high image resolution.

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
Fig. 8 Cross sectional view of a finger and ultrasound images of a finger extensor tendon via conventional signal processing and new signal processing.

Arrows indicate the extensor tendon in each ultrasound image. In the lower image, the new signal processing method filters the signal from the probe to obtain a flat frequency response, as shown in Fig. 7. Here, the high echo lines creating a fibrillar pattern are observed more clearly and in more detail than by conventional signal processing, so that the extensor tendon is more clearly distinguished from its surrounding tissues.

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
Fig. 9 Normal and full screen modes.

A single touch of the monitor display toggles the display between normal and full-screen modes.

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
Fig. 10 Alternative, customizable screen image layouts.

Screen navigation can be customized to a customer's preferences and workflow.

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
Fig. 11 Specialized probes: (a) Hockey-stick linear probe HL18-4 and (b) Micro-convex probe MC10-3.

These probes were specifically developed to access difficult-to-reach parts of the body. The wideband linear probe L11-3 has a view field width of 40 mm. But the HL18-4 probe has a view field width of 25 mm, and the head of the MC10-3 probe has a curvature radius of 14 mm. In spite of such a small head, the MC10-3 probe provides a field angle of 100 degrees or greater.

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
Fig. 12 Comparison between workflows using the X-ray or CT and the ultrasound.

Ultrasound benefits both patients and doctors. Patients needn't move between examination rooms or wait for examination results. Doctors needn't repeatedly open and close medical records, reserve medical equipment, or wait for examination results. Thus, ultrasound improves diagnostic productivity and shortens patients' long and tiring visits to the medical facility.


Fig. 13 Quick and simple retrieval of patient ID.

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