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Emergency Medicine Application of Biomarker Measurement Technology Employing
Surface Plasmon Field-Enhanced Fluorescence Spectroscopy

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Abstract

In the field of emergency medicine, especially in emergency rooms in the United States, the reduction of overcrowding is an urgent task. Patients retained in the ER waste limited hospital resources, and patients suffer from excessively long waiting times for treatment. To reduce overcrowding and provide timely patient treatment, triage must quickly distinguish conditions that require immediate attention, such as acute myocardial infarction (AMI), from less serious conditions whose symptoms (such as chest pain) mimic those of AMI.

In diagnosing AMI, it is recommended to measure the cardiac troponin I (cTnI) protein in the blood with high sensitivity to achieve fast and accurate diagnostic results. At this time, however, such high sensitivity cTnI measurement is possible only in large systems currently found in a hospital's central laboratory, necessitating transport of samples from the emergency room to the central laboratory. Further, since whole blood cannot be used as a measurement sample, the preliminary step of preparing plasma samples from whole blood must be taken. All this takes over 60 minutes for the ER physician to receive results.

To meet this challenge, we developed a small, prototype system that can be kept in the ER, can use whole blood specimens, and can complete high sensitivity cTnI measurement in just 10 minutes by using surface plasmon field-enhanced fluorescence spectroscopy (SPFS).

The analytical sensitivity of this system is equivalent to existing systems in hospital central laboratories that provide high sensitivity cTnI measurement, and agreement of measurements with hospitals' central laboratory systems is good. Furthermore, in a matrix test using whole blood specimens prepared by spiking healthy whole blood with a standard cTnI antigen, a strong correlation between whole blood and plasma was confirmed. The use of this system in the ER promises to greatly aid triage identification of AMI and non-AMI patients, thus reducing ER overcrowding.

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
Fig. 1 High sensitivity cTnI immunoassay system's rapid workflow in diagnosing acute myocardial infarction.

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
Fig. 2 The cTnI SPFS system. The system consists of disposable reagent cartridges and a desktop measuring device.

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Table 1 Specifications of the cTnI SPFS system.

Sample type	Heparinized plasma Heparinized whole blood
Sample volume	200 µl
Turnaround time (TAT)	10 minutes
Measurement range	13 ng/L to 50,000 ng/L (between-run CV < 10%)
Number of measurement slots	3

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Table 2 Analytical sensitivities and within-run imprecisions of our cTnI SPFS system and company X's system for central hospital laboratories.

	Our cTnI SPFS system	Company X's system
LoB	1.1 ng/L	0.8 ng/L
LoD	2.6 ng/L	2.5 ng/L
LoQ	7.2 ng/L	6.9 ng/L
CV 10%		
Within-run imprecision (CV)	52 ng/L 3.8%	N.D.
	315 ng/L 3.8%	N.D.
	48,000 ng/L 2.6%	N.D.


Fig. 3 Dilution test of two cTnI positive heparinized plasma samples with different concentrations spiked with in-house standard cTnI antigen.

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
Fig. 4 Correlation of measurement data: our cTnI SPFS system versus company X's system.

Measurement was conducted on 60 cTnI positive heparinized plasma samples. The range of measured cTnI concentrations was from 3 to 46,173 ng/L.

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
Fig. 5 Correlation of measurement data with our cTnI SPFS system: heparinized plasma samples versus heparinized whole blood samples. Whole blood samples were prepared by spiking in-house standard cTnI antigen to commercially available healthy whole blood. The range of measured cTnI concentration was 4 to 60,256 ng/L.

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