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Ultrasound System SONIMAGE MX1 Realizing Customer Value

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Abstract

In recent years, the market of compact ultrasound diagnostic equipment is expanding in association with improvement of performance and image quality of compact ultrasound diagnostic equipment.

Konica Minolta launched the ultrasound system SONIMAGE HS1 in 2014, and this system has a product specification conforming to the idea of "POC (Point of care)" under the concept of "Visible", "Easy", and "Connectible", and has been accepted in the market.

In 2017, we started to develop a compact and portable ultrasound diagnostic equipment without sacrificing image quality, and launched the ultrasound system SONIMAGE MX1 in March 2018. In the early stage of the development of this system, we collaborated with the design division and the planning division to determine the product concept by turning around a hypothesis-verification cycle. In addition, we repeatedly interviewed target customers about product specifications and designs to reflect them on the product specification.

In this development, interfaces were simplified for improvement of portability and reduction in diagnostic time, so that a user-friendly operation was achieved and convenience in examination was improved.

In addition, image quality and operability have been improved.

Regarding image quality, we have developed a transmission and reception technique, which we call Dual Sonic, and we have achieved improvement of SNR (Signal to Noise Ratio) by about 4 dB due to reduction in an acoustic noise in a shallow area and due to increase in application voltage in association with reduction in unnecessary heat generation and unnecessary acoustic output.

Regarding improvement of operability, we have developed the MPA (Multi Parameter Adjuster) function that optimizes an image, depending on the display depth. Owing to this function, high customizability, and the improved operability, SONIMAGE MX1 won Good Design Award 2018 and receives high valuations from external organizations.

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Fig. 1 Ultrasound system SONIMAGE HS1 available in the market since 2014.

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Fig. 2 Mock-up model of SONIMAGE MX1 at early stage of development.

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Fig. 3 Examples of prototype device made in hypothesis-verification cycles.

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
Fig. 5 Size comparison between SONIMAGE HS1 and SONIMAGE MX1. SONIMAGE MX1 is smaller and lighter than SONIMAGE HS1 by about 40%, and is one hand portable.

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
Fig. 6 Rear view of SONIMAGE MX1 with additional battery.

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Fig. 7 Cradle from which SONIMAGE MX1 can be easily detached when necessary. SONIMAGE MX1 is detached from the cradle by being pulled up in the blue arrow direction.

Fig. 8 SONIMAGE MX1 being used in treatment of knee. SONIMAGE MX1 can be placed behind the patient in line with the affected area.

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
Fig. 9 Basic theories of Triad-THI and Dual Sonic.

(a) Triad-THI: Acoustic noise enters from a shallow part on the outer part of the opening.

(b) Dual Sonic: Triad-THI transmission is performed in the inner part of the opening, and a low frequency transmission is performed in the outer part, so that an acoustic noise in the shallow area can be reduced.

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
Fig. 10 Example of clinical ultrasound image with SONIMAGE MX1.
 (a) Middle finger flexor tendon (left: conventional method right: new method)
 (b) Lumbar spinal spinous part (left: conventional method right: new method)
 Noise in the shallow area is reduced compared with the conventional method, so that clear and fine view is possible.

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
Fig. 11 Conceptual diagram of MPA.

With the MPA function, it is possible to switch to one of image quality settings previously set for each depth by only changing a display depth.

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Fig. 12 Simple Guide on display screen.

Functions and parameters to be set from the console can be checked on Simple Guide displayed on the lower part of the display screen so that operation can be performed without hesitation.

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
Fig. 13 Measurement cursor.

The measurement cursor can be adjusted to be positioned slightly off the touch point so that the cursor position can be finely adjusted without the cursor being blocked by a finger.

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