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High-sensitive Total PSA Test Using Surface Plasmon Field-enhanced Fluorescence Spectroscopy (SPFS) for
Point-of-Care Testing (POCT)

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Abstract

Cases of prostate cancer in Japan have increased rapidly since around 2000. Currently in the diagnosis of prostate cancer, needle biopsies are indicated chiefly by the results of various examinations that determine prostate specific antigen (PSA) values. Since PSA is highly organ specific, it is also effective in predicting the early progression of the disease, making it one of the most widely used tumor markers. A PSA test is widely used as an indicator in prognostic monitoring and recurrence judgment after various treatments. Currently, PSA measurement must be performed with large apparatus, for no compact system has existed that realizes the performance of a large system. Such a compact system would render PSA tests more accessible at clinics and hospitals both in Japan and abroad.

We have created a compact system using surface plasmon field-enhanced fluorescence spectroscopy (SPFS), which measures PSA with high sensitivity by optimizing the reaction protocol for PSA measurement. The system's analytical sensitivity is equivalent to the high sensitivity of existing systems, and measurements made with our SPFS system are in strong agreement with conventional system measurements. In addition, our system provides a wider dynamic range than conventional systems.

Using this new system, clinics whose size and economics previously rendered them unable can now conduct PSA examinations and provide greater individualized medical attention to an increasingly aging society.

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
Fig. 1 PSA testing workflow in diagnosing prostate cancer.

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
Fig. 2 Compact high-sensitivity PSA measurement system consisting of disposable reagent cartridges, disposable contamination-free resin pipette tips, and a desktop measuring device.

The PSA reagent cartridge is mounted with an SPFS optical prism chip on which an anti-PSA antibody is solidified, and is filled with other reagents necessary for measurement (i.e. sample dilution liquid, fluorescent labeled antibody liquid, and cleaning liquid). The operator attaches the measurement reagent cartridge to the desktop measuring device and dispenses the sample into the cartridge to start measurement.

Table 1 PSA SPFS system specifications.

sample type	serum
sample volume	200 μ L
Turnaround time (TAT)	10 minutes
Number of measurement slots	3
Temperature	15 – 30 $^{\circ}$ C

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Table 2 Reagent handling performance of PSA SPFS system.

The PSA SPFS system achieves highly precise measurement reproducibility with highly precise, automated reagent handling for the dilution of samples.

Dilution rate	6 dilution \pm 0.22 dilution
Dispensing accuracy	150 μ L \pm 6 μ L
Dispensing repeatability (CV)	1.2 %
Residual liquid volume	< 7 μ L

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Table 3 Analytical sensitivities and within-run imprecisions.

The PSA SPFS system exhibits LoB (limit of blank), LoD (limit of detection), and LoQ (limit of quantification) equivalent to those of the PSA CLEIA method, as well as low within-run imprecision CV (coefficient of value) and wide measurement dynamic range.

	Our PSA SPFS system	PSA CLEIA method
LoB	0.0023 ng/mL	0.0016 ng/mL
LoD	0.0054 ng/mL	0.0030 ng/mL
LoQ (20%CV)	0.0055 ng/mL	0.0043 ng/mL
Within-run imprecision	0.1 ng/mL 5.1% CV	N.D.
	1.0 ng/mL 4.6% CV	N.D.
	70 ng/mL 5.2% CV	N.D.
Measurement range	0.008–5,000 ng/mL	0.008–100 ng/mL

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
Fig. 3 Calibration curves of total PSA in a standard antigens sample from the SPFS system (●) and company X's CLEIA system (△). The SPFS system has a wider measurement range.

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
Fig. 4 Correlative evaluation of SPFS system and company X's system. Measurements were conducted on 48 serum samples whose concentrations were lower than 100 ng/mL total PSA.


Fig. 5 Correlative evaluation of SPFS system and company X's system. Measurements were conducted on 54 serum samples, including 6 cases whose total PSA concentration was over 100 ng/mL.

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Table 4 Contents of six automated immunoassay systems that perform total PSA tests.

Instrument	Company	Size (mm) (W)×(D)×(H)	Measurement range (ng/mL)	TAT
-	Konicaminolta	450×490×490	0.003-5000	10 min
ARCHITECT	Abbott	1549×1245×1219	0.008-100	29 min
UniCel	Beckman	1710×970×1700	0.008-150	20 min
cobas	Roche	1200×730×800	0.002-100	>18 min
AIA	Tosoh	1960×840×1360	0.003-100	18 min
POCube®	Toyobo	280×310×275	0.1-80	6 min

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
Fig. 6 Rapid workflow of a high-sensitivity PSA immunoassay system in diagnosing prostate cancer.