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Development of Patient Dose Management System FINO.XManage

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Abstract

Due to the record-keeping and management of medical exposure dose becoming mandatory in April 2020, exposure dose management systems are increasingly demanded. We have leveraged our strengths as a manufacturer of imaging diagnostic-related equipment, such as PACS (Picture Archiving and Communication System) and computed radiography, and have developed the patient dose management system FINO.XManage, which can manage medical exposure dose, control exposure dose, and improve efficiency of exposure dose management work by linking exposure dose with radiographic images.

Due to labor shortage in the medical site and need for cost reduction, radiation facilities face the need for further improvement in efficiency of their work. In order to carry out exposure dose management appropriately without changing the current work procedure, we have improved convenience in operation in facilities by covering different imaging devices made by different manufactures and by managing a wide range of examination from CT to nuclear medicine tests in an integrated manner. Also, we have mounted some functions strongly considering exposure dose reduction without image quality loss so as to contribute to quality improvement of medical service.

When making product specifications, we have repeatedly performed a cycle of prototype/hearing/feedback with many specialists (doctors and radiological technologists) from the design stage.

The following technologies will be also described here.

1. Technology for the automatic calculation function of SSDE (Size-Specific Dose Estimates)
2. Function for linking with our image quality assurance system I-PACS QA

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
Fig. 1 Product Concepts of FINO.XManage.

The technologies are established to achieve the following three production concepts.

1. Openness of the system such that the system can be proposed for small to large facilities.
2. Achievement of both exposure dose reduction and no image quality loss.
3. Improvement in efficiency of work by automatic determination processing.

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
Fig. 2 Configuration diagram of FINO.XManage System

The basic configuration includes a DICOM communication section, a database, and a web server.

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Table 1 Integrated processing of exposure dose data.
In medical sites, there are many conventional medical devices which RDSR is not available or incomplete. The FINO.XManage is equipped with a process to integrate RDSR data so that all types of exposure dose data can be managed in an integral manner.

Case classification	Example functions of FINO.XManage
RDSR not available	Exposure dose data is automatically generated by performing OCR of a dose summary image (e.g. a captured image of a monitor display) of CT.
	RDSR is automatically generated from DICOM header information of a mammography image.
	Nuclear medicine dose information is acquired in CSV format from RIS.
	Exposure dose data of an angiography (XA) is generated with user's manual input.
RDSR information incomplete	Height and bodyweight information is automatically acquired, for supplement, from a DICOM image or an information system.
	Protocol data is obtained by automatically analyzing information in a DICOM header and is set to RDSR.

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
Fig. 3 Scatter plot (CTDIvol-bodyweight).
Horizontal axis: weight
Vertical axis: CTDIvol
High operability is achieved such that a high-quality graph can be generated by one click.

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
Fig. 4 Flowchart of examination with high exposure.
A technique for analyzing exposure dose for each imaging protocol, based on patient attributes such as age and bodyweight, and it is made easy to switch from a scatter plot of exposure dose or the like to exposure dose information for the examination, thereby improving convenience in an analysis work of an examination with high exposure.

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Table 2 Image retention period depending on purpose.

To reduce running cost, on the side of facilities, for storing medical images, images are compressed or thinned depending on the purpose of use, and the retention period can be individually set.

Purpose	Retention Period
JPEG files for Viewer	as far as possible
DICOM files for SSDE, etc.	2 weeks
DICOM files for SD value	6 months

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
Fig. 5 Operations in cooperation with I-PACS QA.

Konica Minolta image quality assurance system I-PACS QA can display, on the display screen of an examination image, exposure dose information and an alert in case of overexposure.

Fig. 6 Alert marks in a CT examination in I-PACS QA.

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
Fig. 7 Calculation of SSDE (According to AAPM Report No. 204).

In CT examinations, there is used SSDE (Size-Specific Dose Estimates) that takes into consideration the difference in physique between patients. For FINO.XManage, there has been developed a function for calculating SSDE by calculating the Lateral (LAT) dimension and the AP dimension of an abdominal region from an Axial image of a CT examination.

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- 1) AAPM Report No.204 (2011) : Size-Specific Dose Estimates (SSDE) in Pediatric and Adult body CT Examinations