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Highly Sensitive Quantitative Evaluation of HER2 by Ultrasensitive Tissue Staining Service (Quanticell®) Using Fluorescent Nanoparticle PID (Phosphor Integrated Dots)

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Abstract

In order to promote personalized medicine, there is a need to develop an *in vitro* diagnostic drug as a criterion for selecting an appropriate treatment method. Thus, in recent years, development of companion diagnostic drugs (CDx) linked to specific drugs has been accelerating.

Human epidermal growth factor receptor 2 (HER2) is over-expressed in about 15-20% of metastatic breast cancers and is a predictor of the efficacy of the anti-HER2 therapeutic agent trastuzumab. The immunohistochemistry (IHC) HER2 protein detection method has been approved as a companion diagnostic (CDx) for anti-HER2 therapeutics. However, since the approved CDx is pigmentation-based DAB (3,3'-diaminobenzidine) stain, the judgment is semi-quantitative by visual inspection. Furthermore, since the signal is amplified by the peroxidase enzyme reaction, it is difficult to judge when the expression level is low.

Low-level HER2 expression cases without HER2 amplification are observed in about 50% of breast cancers. However, recently developed anti-HER2 ADC therapies have shown tumor shrinkage even in patients with low HER2 expression (HER2 IHC 1+). Therefore, there is no diagnostic method for highly accurate stratification of patients with low HER2 expression.

To solve this problem, by using the Quanticell® ultra-sensitive tissue staining service with PID (Phosphor Integrated Dots) which are high-intensity fluorescent nanoparticles originally developed by Konica Minolta, we introduce an example of detecting HER2 protein in formalin-fixed paraffin embedding (FFPE) tissues.

It is shown that Whole Slide Image (WSI) with Quanticell® displays HER2 protein even when the clinical tissues are classified as HER2-negative by the current IHC method. In addition, it is shown that stratification of HER2-low-expressing patients from HER2-negative patients is possible.

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Table 1 The FDA-approved HER2 IHC test kits HercepTest Dako, Oracle Leica, and PATHWAY Ventana are used to properly evaluate anti-HER2 drug treatment cases.

All existing test kits are based on the DAB method and are suitable for stratifying HER2 3+/2+ with high HER2 expression levels. However, they are unsuitable for stratifying HER2 1+/0.

Drug	Antigen	1stAb	Type	Company	Detection	Frequency
Trastuzumab	HER2/Neu	A0485	Polyclonal	DAKO	DAB	21%
Trastuzumab	HER2/Neu	4B5	Monoclonal	Ventana	DAB	63%
Trastuzumab	HER2/Neu	CB11	Monoclonal	Leica	DAB	5%

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
Fig. 1 Schematic diagram of the Quanticell[®]-based assay and comparison of DAB and PID staining with TMA.

The particle size of the fluorescent nanoparticles PID is about 130 nm (excitation wavelength: 580 nm, emission wavelength: 620 nm). Using 104 cases / 208 cores (US Biomax, BR20810) of breast tumor tissue FFPE microarray (TMA), staining was performed using DAB stain (Ventana I-VIEW), as well as the same primary antibody (clone 4B5) as DAB stain. Staining was also performed using Konica Minolta Fluorescent Nanoparticles (PID) in place of DAB with secondary antibodies.

TMA scoring was performed according to the ASCO / CAP HER2 test guidelines, and each core was classified into 3+/2+/1+/0 based on DAB staining. The TMA results for DAB staining and PID staining were in good agreement overall.

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
Fig. 2 Assay based on Quanticell®, using Konica Minolta's new fluorescent immunohistochemical technology.

Since PID staining has higher sensitivity than DAB staining, it can clearly detect HER2 proteins with extremely low expression levels. In addition, because it has a wide dynamic range, both HER2 over-expressing (3+) and low expressing (1+/0) samples are identified under the same staining conditions even with TMA containing many tissues with different HER2 expression levels.

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
Fig. 3 Comparison between DAB staining and PID staining in tissues with different HER2 IHC scores (3+, 2+, 1+, 0).

Since DAB staining of HER2 1+/0 shows a weak signal, it is not easy to make an objective and reproducible judgment. However, with PID staining, a clear positive signal can be detected even in cores classified as HER2 1+ by the DAB method. Circumferential HER2 membrane staining is observed with cores which are classified as HER2 0 expression with DAB staining.

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

Fig. 4 A PID-stained image of tissues judged as low HER2 expression by DAB staining.

The PID value is an integral of the fluorescence intensity of PID in the invasive cancer area. Clear positive signals were detected in all cores, classified as HER2 1+, which is difficult to determine by DAB staining.

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Fig. 5 Comparison between DAB and PID scores. The PID-stained image was converted into a high-precision Whole Slide Image (WSI) by a virtual slide scanner, and the integrated fluorescence intensity was extracted from the invasive cancer area.

In accordance with ASCO / CAP HER2 test guidelines, the average fluorescence intensity of the top 10% was calculated as the PID score. Since Quanticell® calculates the score based on the fluorescence intensity, a linear score can be obtained, which is continuous rather than a discontinuous score such as 3+ / 2+ / 1+ / 0.

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■ = 0+ by DAB IHC but considered 1+ by Quanticell
 ■ = 0+ by DAB IHC and Quanticell

Table 2 Comparison between DAB classification and estimated Quanticell® classification.

A concordance table of PID scores and DAB scores was created. In the HER2 3+/2+ category, the Quanticell® and DAB methods were in very good agreement. On the other hand, due to the high-sensitivity detection of Quanticell®, about 40% of the HER2 1+ cores show a positive signal corresponding to HER2 2+. Also, for the HER2 0 core, a positive signal corresponding to HER2 1+ was observed at a rate of around 50%.

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