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Introduction of Intelligent Workstation “KINOSIS” for Dynamic Digital Radiography

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Abstract

In order to overcome the limitation on the diagnostic capability of chest X-ray image, Konica Minolta has developed new imaging modality named Dynamic Digital Radiography (DDR). We make DDR ready with simple configuration like a conventional chest X-ray images, which is a static image and widely used as a first-line imaging test in primary care. Compared to conventional chest X-ray images, DDR provide a dynamic image, thus, it can increase the amount of information in initial diagnosis phase. We launched “KINOSIS” (English name: “Intelligent Workstation”) in 2018 which is a dedicated imaging workstation for DDR.

Our Intelligent Workstation provides medical fields with high accuracy of diagnosis by “visualization of movement” for better visibility and “quantification of movement” by digitizing and displaying graphs of the movement. It is utilizing our image processing technology which was developed for conventional radiography analysis.

After the launch, gathering opinions about time constraints at medical practices and differences of needs among different clinical departments, we have improved the system so that it can access the information quickly and multidirectionally. Furthermore, we are introducing new functions, such as information integration to the system, and working on enhancing accuracy and expanding the clinical applications by collaborating with domestic and overseas hospitals.

This paper reports the features of the dynamic digital radiography system and examples of clinical applications, created from simple X-ray imaging, together with the latest trends in the market.

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
Fig. 1 Overview of hospital IT system.

Intelligent Workstation is an X-ray imaging workstation in hospital IT system, and it manages analysis and output of dynamic X-ray image.

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
Fig. 2 Software structure of automatic analysis processing of Intelligent Workstation.

Intelligent Workstation can fully automatically execute all the multifaceted image analysis processing, required for multiple clinical departments by only one examination.

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
Fig. 3 Software structure to control image distribution route of Intelligent Workstation.

Intelligent Workstation can distribute processed outputs to not only PACS but entire hospital IT system as web contents.

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
Fig. 4 Schematic view of variation of image analysis on Intelligent Workstation.

Intelligent Workstation provides realistic and easy observation by frequency processing and structure recognition processing on dynamic X-ray image. By using image analysis technology, it is enabled to visualize and quantify respiratory and circulatory organ-related biological and physiological functions.

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
Fig. 5 PH2-Mode algorithm.

In PH2-Mode, end-diastole frame is used as reference, and the amount of change in the signal, synchronized with the heartbeat waveform, is obtained by taking the difference between each frame and the reference. Quantitative evaluation of blood circulation is expected.

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
Fig. 6 Comparison of visibilities for pulmonary perfusion defects.

Left: PH-MODE, middle: PH2-MODE, right: pulmonary perfusion scintigraphy.

Compared to the PH-MODE, PH2-MODE has improved visibility for pulmonary perfusion defects. In addition, an image similar with pulmonary perfusion scintigraphy is obtained.

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Fig. 7 Intelligent Workstation: structure observation mode.

This mode supports surgical planning by offering information about pleural adhesion and invasion. The abnormality of tissue motion, which suggests adhesion and invasion, is visualized with FE-MODE and DM-MODE.

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Fig. 8 Intelligent Workstation: dynamic image evaluation mode.

This mode supports treatment planning/evaluation by offering information about quantified motions of lung and diaphragm mainly to respiratory department. Automatically measured diaphragm excursion and lung field area at specific timing as well as a trend chart compared to past exams are used.

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- 1) Rie Tanaka, Tohru Tani, Norihisa Nitta, et al., "Detection of Pulmonary Embolism Based on Reduced Changes in Radiographic Lung Density During Cardiac Beating Using Dynamic Flat-panel Detector: An Animal-based", Study. Academic Radiology, 26 (10), 1301 (2019)
- 2) Rie Tanaka, Shigeru Sanada, Masaki Fujimura, et al., "Development of pulmonary flow evaluation method with a dynamic flat-panel detector: quantitative correlation analysis with findings on perfusion scan", Radiological Physics and Technology, 3, 40 (2010)
- 3) 松谷 哲嗣, "JIRAテクノロジーレポート", 57, 30, (2018)

- 4) Yuzo Yamasaki, Kohtaro Abe, Kazuya Hosokawa, Takeshi Kamitani, “A novel pulmonary circulation imaging using dynamic digital radiography for chronic thromboembolic pulmonary hypertension”, *European Heart Journal*, 41 (26), 2506 (2020)
- 5) 高田 宗尚, 第2回X線動態画像セミナー第2部臨床研究報告-X線動態解析による肺血流評価: 肺血流シンチグラフィとの比較, *INNERVISION*, 35 (3), 2 (2020)
- 6) X線動画解析ワークステーションKINOSIS—紹介動画・ドクターズインタビュー—<<https://www.konicaminolta.jp/healthcare/products/dr/kinosis/index.html>>
- 7) N. ohkura, et. al., “Dynamic-ventilatory digital radiography in air flow limitation: change in lung area reflects air trapping, *Respiration*”, 99, 382, (2020)
- 8) T. Hida, et. al., “Decreased and slower diaphragmatic motion during forced breathing in severe COPD patients: Time-resolved quantitative analysis using dynamic chest radiography with a flat panel detector system, *European Journal of Radiology*”, 112, 28 (2019)
- 9) 大倉 徳幸, 第2回X線動態画像セミナー第2部臨床研究報告—X線動態解析を用いた新たな呼吸機能評価の検討, *INNERVISION*, 35 (3), 5 (2020)
- 10) T. Hino, et. al., “Projected lung areas using dynamic X-ray (DXR), *European Journal of Radiology Open*”, 7, 100263 (2020)