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Development of Ultrasound Diagnostic System SONOVISTA GX30 for Obstetrics and Gynecology

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Abstract

Konica Minolta has placed on the market ultrasound diagnostic system suited for POC (Point of care, diagnosis on the site) under the concept of "Visible", "Easy", and "Connectible", and the system has been highly evaluated in the field of orthopedics for its high-definition image and intuitive operation.

Meanwhile, also in the field of obstetrics and gynecology, we are expanding the business field by providing ultrasound diagnostic system suited to the workflow of obstetricians and gynecologists under the same concept.

We have developed this time the SONOVISTA GX30, which has a compact body and flexibility to fit in various layouts and can be installed flexibly even in a small examination room.

This paper describes various technologies for realizing this ultrasound diagnostic system SONOVISTA GX30, which is valuable for obstetricians and gynecologists and pregnant women.

Regarding the concept of "Visible", we have established a technology, based on a new virtual origin method that provides a wide viewing angle and a high image quality in the center part of the image that obstetricians and gynecologists watches most carefully. In addition, we have independently developed an image processing technology to convert a low definition image into a high definition image while keeping artifacts low, so that the diagnostic performance has been improved.

Regarding the concept of "Easy", We have developed technologies to achieve improvement in operability of a console and a touch panel, easy measurement of fetal growth, and a function for directly storing fetal images on an external recording medium.

As a result, Konica Minolta's identity has been combined with easy operability of currently used ultrasound diagnostic system to achieve easy-to-use functions and easy operation in the field of obstetrics and gynecology from the point of view of measurement and various operations.

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Fig. 1 Obstetrical and gynecological examination room and internal examination room.

Because various kinds of instruments necessary for diagnosis are placed in an internal examination room, an ultrasound diagnostic system needs to have a small footprint for flexible installation.

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Table 1 Footprint of the SONOVISTA GX30.

The housing is made compact, and the center of gravity is lowered, so that the footprint is made equal to or smaller than that of the current models.

	Footprint (cm ²)	Width-depth range (mm)
SONOVISTA GX30	2,065	W 386 x D 535
SONOVISTA FX PE (Current equipment)	2,112	W 389 x D 543

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
Fig. 2 Movable range of the SONOVISTA GX30.

The console of the SONOVISTA GX30 can be rotated 90 degrees leftward and 180 degrees rightward, and the console monitor can be moved vertically by 120 mm.

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Table 2 Movable range.

We have achieved such a structure that the SONOVISTA GX30 has a larger movable range than the current models such that the SONOVISTA CX30 is more optimally positioned vertically and rotationally than before depending on obstetrician and gynecologist use.

	Swivel range (deg)	Up-down range (mm)
SONOVISTA GX30	270	120
SONOVISTA FX PE (Current equipment)	90	100

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Fig. 3 Endovaginal probe to be used with the SONOVISTA GX30.


Fig. 4 Effect of a wide viewing angle of 220 degrees, which is a development target.

The right and left ovaries can be simultaneously diagnosed with a wide viewing angle of 220 degrees, which is the development goal.

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
Fig. 5 Difference in B-mode image quality of a fetal brain between with and without iXRET (eXtended Resolution Enhancement Technology). Fig. 5(b) shows a fetal cranium more clearly with iXRET applied than Fig. 5(a).


Fig. 6 Difference in quality of color flow mode images between with and without iXRET. Fig. 6(b) shows a blood flow in more detail with iXRET than Fig. 6(a).

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Fig. 7 Operation on console and touch panel of SONOVISTA GX30.

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
Fig. 8 Two measurement modes of SONOVISTA GX30.

- (a) A fetal sac length (distance) is first measured, and a fetal sac label is then selected.
- (b) A fetal sac label is first selected, and a fetal length is then measured.

The both modes are available with SONOVISTA GX30 to cover many obstetricians' and gynecologists' ways of use.

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